


Graphite Tubes for AAS



I. Product Line:

Our graphite tubes for AAS include six models as follows:

- Longitudinal-heating standard graphite tube
- Longitudinal-heating pyrolytically-coated graphite tube
- Longitudinal-heating pyrolytically-coated graphite tube with platform
- Transversal-heating pyrolytically-coated graphite tube with platform
- Transversal-heating high-sensitivity graphite tube
- Transversal-heating high-sensitivity graphite tube for Jena


II. Design:

Our products are specially designed to meet or exceed the requirements of various types of atomic absorption spectrophotometers, providing long lifetime for your most difficult analysis. We customize and OEM all types of graphite tubes.


III. Quality Control:

The components are made from high-purity graphite, and each one is carefully examined prior to shipment. We guarantee 100% customer satisfaction, high quality and low prices.


IV. List of Graphite Furnace Products:



Type: Longitudinal-heating standard graphite tube for Perkin Elmer AAS



Type: Longitudinal-heating coated graphite tube for Perkin Elmer AAS




Type: Longitudinal-heating coated graphite tube with platform for Perkin Elmer AAS




Type: Transversal-heating coated graphite tube with platform for Perkin Elmer AAS




Type: Transversal-heating coated high-sensitivity graphite tube with platform for Perkin Elmer AAS




Type: Longitudinal-heating standard graphite tube for Hitachi AAS




Type: Longitudinal-heating coated graphite tube for Hitachi AAS




Type: Longitudinal-heating coated graphite tube with platform for Hitachi AAS




Type: Longitudinal-heating standard graphite tube for Shimadzu AAS




Type: Longitudinal-heating coated graphite tube for Shimadzu AAS



Type: Longitudinal-heating coated graphite tube with platform for Shimadzu AAS




Type: Longitudinal-heating standard graphite tube for Varian AAS




Type: Longitudinal-heating coated graphite tube for Varian AAS




Type: Longitudinal-heating coated graphite tube with platform for Varian AAS




Type: Longitudinal-heating standard graphite tube for Unicam AAS




Type: Longitudinal-heating coated graphite tube for Unicam AAS



Type: Longitudinal-heating coated graphite tube with platform for Unicam AAS




Type: Longitudinal-heating standard graphite tube for GBC AAS




Type: Longitudinal-heating coated graphite tube for GBC AAS



Type: Longitudinal-heating coated graphite tube with platform for GBC AAS



Type: Transversal-heating coated graphite tube with platform for GBC AAS



Type: Transversal-heating coated high-sensitivity graphite tube with platform for JENA NOVA AAS (deuterium lamp background correction)



Type: Transversal-heating coated high-sensitivity graphite tube with platform for JENA Zeenith AAS (Zeeman background correction)



Type: Longitudinal-heating coated graphite tube for Thermo AAS



Type: Transversal-heating coated graphite tube with platform for Purkinje AAS



Type: Transversal-heating coated graphite tube with high-sensitivity for Purkinje AAS



Graphite Electrode

For many manufacturers, like GBC, Perkin Elemer, Shimadzu, Unicam, Varian and etc.